

TECHNICAL SPECIFICATIONS

 OCTOBridge™	Combined Pack			Implant Only		
	1.4 OCTOBridge™	1.7 OCTOBridge™	2.4 OCTOBridge™	1.4 OCTOBridge™	1.7 OCTOBridge™	2.4 OCTOBridge™
SKU	OCBR14N	OCBR17ND	OCBR24ND	OCBR14Ni	OCBR17Ni	OCBR24Ni
Description	OCTOBridge Push 1.4 with Needles	OCTOBridge Push 1.7 with Needles and Cannulated Drill	OCTOBridge Push 2.4 with Needles and Cannulated Drill	OCTOBridge 1.4 with Needles	OCTOBridge 1.7 with Needles	OCTOBridge 2.4 with Needles
K-wire size	1.4mm	1.0mm	1.2mm	1.4mm	1.0mm	1.2mm
Drill size	1.4mm	1.75mm	2.4mm	1.4mm	1.75mm	2.4mm
Drill Depth	12mm	12mm	17mm	12mm	12mm	17mm
Suture USP	#2-0	#0	#2	#2-0	#0	#2
No. of Needles	2	2	2	2	2	2
Needle type	15mm, 3/8 circle	22mm, 3/8 circle	22mm, 3/8 circle	15mm, 3/8 circle	22mm, 3/8 circle	22mm, 3/8 circle
 OCTOExtend™	Combined Pack			Implant Only		
		1.7 OCTOExtend™	2.4 OCTOExtend™		1.7 OCTOExtend™	2.4 OCTOExtend™
SKU	-	OCEX17D	OCEX24D	-	OCEX17i	OCEX24i
Description	-	OCTOExtend Push 1.7 with Cannulated Drill	OCTOExtend Push 2.4 with Cannulated Drill	-	OCTOExtend 1.7	OCTOExtend 2.4 and Suture Loop with Needle
K-wire size	-	1.0mm	1.2mm	-	1.0mm	1.2mm
Drill size	-	1.75mm	2.4mm	-	1.75mm	2.4mm
Drill Depth	-	12mm	17mm	-	12mm	17mm
Suture USP	-	#2-0	#2-0	-	#2-0	#2-0
No. of Needles	-	N/A	N/A	-	N/A	N/A
Needle type	-			-		
 OCTOAnchor™	Combined Pack			Implant Only		
	1.4 OCTOAnchor™	1.7 OCTOAnchor™	2.4 OCTOAnchor™	1.4 OCTOAnchor™	1.7 OCTOAnchor™	2.4 OCTOAnchor™
SKU	OCAN142N	OCAN172ND	OCAN242ND	OCAN142Ni	OCAN172Ni	OCAN242Ni
Description	OCTOAnchor Push 1.4 Double Loaded with Needles	OCTOAnchor Push 1.7 Double Loaded with Needles and Cannulated Drill	OCTOAnchor Push 2.4 Double Loaded with Needles and Cannulated Drill	OCTOAnchor 1.4 Double with Needles	OCTOAnchor 1.7 Double with Needles	OCTOAnchor 2.4 Double with Needles
K-wire size	1.4mm	1.0mm	1.2mm	1.4mm	1.0mm	1.2mm
Drill size	1.4mm	1.75mm	2.4mm	1.4mm	1.75mm	2.4mm
Drill Depth	12mm	12mm	17mm	12mm	12mm	17mm
Suture USP	#2-0	#0	#2	#2-0	#0	#2
No. of Needles	4	4	4	4	4	4
Needle type	15mm, 3/8 circle	22mm, 3/8 circle	22mm, 3/8 circle	15mm, 3/8 circle	22mm, 3/8 circle	22mm, 3/8 circle

