

THE NX NAIL TECHNICAL SPECIFICATIONS



Nail Size	Indication	Poximal Head \varnothing	Shaft \varnothing	Distal Tip \varnothing	Shortest Implant	Longest Implant	K-Wire \varnothing	Drill \varnothing	Screw Driver
2.0	Phalanx	2.7	1.5	2.0	12	40	0.8	1.5	1.5mm Hex
2.5	Phalanx	3.15	1.75	2.5	12	60	1.0	1.75	1.7mm Hex
3.0	Metacarpal	3.9	2.0	3.0	25	70	1.0	2.0	T8 Torx
3.5	Metacarpal	4.8	2.4	3.5	25	70	1.2	2.4	T10 Torx
4.0	Metacarpal	5.5	2.75	4.0	35	80	1.2	2.75	T15 Torx
5.0	Metacarpal	6.9	3.5	5.0	35	90	1.6	3.5	T20 Torx

Note: Measurements are shown in mm

